


LokTab™ Ordering Guide

The LokTab™ Patient Notifier can be ordered in three distinct configurations with additional options per configuration. Use the below guidelines to build the notifier to best meet your unique needs.

LokTab Type A

The header panel is adhered to the face of LokTab™ Notifier.
The only collection available for this option is Phoenix (PH).


- 1 **How many sliding tabs are needed?**
 - 8 Tabs (4 Left / 4 Right)
 - 10 Tabs (5 Left / 5 Right)
- 1 **Which configuration would you like for the face?**
 - Full Window
 - 2 Half Windows
 - Full Whiteboard
 - 1 Window / Whiteboard
 - Full Accent
- 3 **Would you like to add an accessory?**
 - Aluminum Note Holder
 - Decorative Footer (Phoenix)



LokTab Type B

The header panel is mounted directly above LokTab™ Notifier,
with a small gap between the two.

- 1 **How many sliding tabs are needed?**
 - 8 Tabs (4 Left / 4 Right)
 - 10 Tabs (5 Left / 5 Right)
- 2 **Which collection would you like for the Header?**
 - Aspen
 - Rio
 - Phoenix
 - Sydney
 - Florence
- 3 **Which configuration would you like for the face?**
 - Full Window
 - 2 Half Windows
 - Full Whiteboard
 - 1 Window / Whiteboard
 - Full Accent
- 4 **Would you like to add an accessory?**
 - Aluminum Note Holder
 - Decorative Footer (Aspen or Phoenix)



LokTab Type C

The full face of the LokTab™ Notifier is covered by the sign panel.
The only collection available for this option is Phoenix (PH).

- 1 **How many sliding tabs are needed?**
 - 8 Tabs (4 Left / 4 Right)
 - 10 Tabs (5 Left / 5 Right)
- 2 **Would you like to add an accessory?**
 - Aluminum Note Holder
 - Decorative Footer (Phoenix)

